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STEP 1

TURN ON THE DEVICE

(the ON/OFF button is located on the top of the back cover).

SCAN NFC CODE

(using the NFC code you will gain access to the electronic manual (pdf).

Link to the manual:

<https://www.ensoel.pl/support/sweeper-ct/>

STEP 2



MAKE THE DECISION to choose a program with:

A. Quick Access Buttons (requires programming)

See page 14

B. MY PROGRAMS (requires programming - manually from the Base Catalog or the optionally attached professional program LIFE COMMANDER 4) - Patrz strona 14

C. BASE CATALOG (ready-made perfectly matched 663 programs)

STEP 3

SELECT A READY PROGRAM from BASE CATALOG:


1. The highlighted/selected program changes the color of the letters to yellow.
2. To make sure how many programs have been selected, look at the numbers on the right side of the RUN button.
3. Press the RUN button.

STEP 4

ELECTRODE CONNECTION TEST

CHECK IF THE ELECTRODES ARE PROPERLY CONNECTED.

The white (previously red) arrows indicate where you need to connect the cables (5) and electrodes (6) in a given program. After correctly connecting the electrodes, the color of the arrows will change to green. Now you can press the START button.

STEP 5

CLICK START IN THE WORK WINDOW:

If you have decided that you do not want to use a particular channel or electrode, the program will ask if you want to disable the channel (electrode) so that the device does not monitor circuit continuity during program execution.

STEP 5a

SCREEN SAVER

To ensure that the battery lasts as long as possible, the screen saver is activated after a few minutes of program use. To reactivate the screen, simply tap the screen once with your finger. In "pocket mode", it takes three taps.

STEP 6

TURNING OFF THE DEVICE

After executing the program, the device will turn off automatically after a few minutes. If you want to turn off the device earlier or exit the WORK WINDOW while the program is running, pause the program by pressing the PAUSE button and then the EXIT icon.

FREQUENTLY ASKED QUESTIONS

TROUBLESHOOTING

Remember: The device is very user-friendly, but very complex in its construction.

1. Why is there a pause during program execution and the device prompts for program continuation?

- Each user will subjectively feel impulses on their palm skin while using the device with electrodes. The device is equipped with default automatic conductivity control (it can be disabled by the "Gel Electrodes" mode in the settings if needed). If at any point the conductivity decreases or there is insufficient moisture or insufficient contact with the skin, the program will pause its operation. A message saying "Electrodes disconnected" will be heard. The possible causes for this could be:
 - a) You should check if the sockets, cables are damaged or if the cables are properly connected to the device.
 - b) Too dry, damaged hand skin or thick calluses.

how to deal with it?

- 1) Lightly dampen your hands and repeat the program.
- 2) Increase the amplitude for the electrodes and repeat the program.
- 3) Activate the "Gel Electrode" mode located in SETTINGS.

REMEMBER: do not use electrodes where the skin in the area of contact with the electrodes is damaged.

2. How to use a loop / induction mat?

- The induction loop affects the entire body regardless of its location, however, to intensify the effects, it can be placed in a diseased or painful area. It can be placed around the neck, draped over the shoulder, placed on the hips, joints, etc. (when connecting an external induction loop, the internal induction automatically turns off).
- The induction mat is used either as a backrest or can be sat upon on a chair. Due to its design, the induction mat provides a more uniform field, ensuring better parameters than with the induction loop.
- The internal induction always works when any program is activated on the device, so frequency sessions can be used without connecting any electrodes or external induction loops. Just set the desired program and carry the device with you. The device can be applied to the area of pain emission or where discomfort is present.

3. What frequencies are covered by the available programs?

- Here, a wide range of possibilities derived from the scientific research of Rife, Clark, as well as other scientists and our own studies, has been applied.

4. Will the general program work for everyone?

- In the general program, there is a wide range of frequencies that have a positive effect in most cases.

5. Are foot electrodes necessary in a sweeper?

- To enable the flow of electricity throughout the entire body and enhance the effects of the session, it is possible to connect foot electrodes. In order to do this, it is necessary to place an individual order for them with Enso Electronics company.

6. What can be the cause of inefficiency?

- If we have been using a device for some time and the results are not satisfactory, it may be related to maintaining an improper diet, which is an essential element during Sweeper sessions, or incorrectly selecting the programs by ourselves. In such a situation, we recommend visiting a specialist/naturopath who uses the Volla or Vega test method, who will be able to perform an electroacupuncture test and tailor the therapy to individual needs.

7. How many sessions per day can be used with a sweeper?

- The number of sessions depends primarily on individual well-being and vitality.
- It is necessary to use supportive measures to cleanse the body and remove toxins generated during therapy.
- Regardless of feeling good, it is not advisable to exceed 4 sessions per day targeting pathogens. This does not apply to programs such as detox, drainage, organ support programs, etc.
- It is recommended to take breaks between sessions. It is recommended that a session does not last longer than 2 hours.

8. What is the order of conducting health therapies?

- The order of actions should be determined by a specialist or therapist. In the absence of such consultation, programs should be set based on symptoms, such as headache, swelling, bloating, etc.

9. Does the device work similarly to a zapper?

- The part generating current is based on the principles of Hulda Clark's zapper, but it has an improved mechanism that primarily utilizes selective frequency ranges, which do not harm the natural microflora and reach the body more effectively.

10. Can the device be used on children?

- It is not recommended to independently use programs on children under the age of 3. It is always advisable to consult with a bioresonance therapist. For children, programs are used with a time of 1 program shortened to 7 minutes.

11. Sensitivity to electric current (feeling electric current)

Many factors influence sensitivity to electric current:

1. Firstly, it depends on the **physical parameters of the signal on the electrodes**:

- It depends on the frequency of the current flow direction changes. The lower the frequency (less frequent current jumps), the greater the sensitivity. One program may consist of multiple signals with different frequencies.
- It depends on the shape of the signal - the steeper the signal, the more it is perceived.
- Of course, it depends on the intensity - the higher, the more noticeable.
- It depends on the polarization of the applied signal.

2. Secondly, the sensitivity is influenced by the **conductivity of the human body** part located between the electrodes. Conductivity also depends on many factors.:

- skin moisture,
- general health condition of people,
- mental state of a person,
- blood pressure,
- weather conditions,
- procedure time,
- number of procedures per day,
- age,
- gender
- and many others.

One cannot occasionally compare the sensation, as external factors often change quite frequently. At different times, a person may feel an electric current or not feel it at all. This is not a sign that the device is not working. It is a normal process based on biophysical laws. The device is designed in such a way that it has controlled electrical parameters that do not change during its operation. In case of any doubts, it can always be checked at the manufacturer's service.

Why is it better not to feel tingling? Familiarize yourself with the Arndt-Schultze rule below.

The Arndt-Schultz law states that a weak stimulus produces a beneficial reaction in the organism, whereas a strong stimulus can negate this effect. It has been shown that relatively low intensities of magnetic fields are more effective than higher intensities.

The use of significantly higher intensity fields (based on the principle of "more is better") is undoubtedly ineffective and may even have the opposite effect. It is also important to consider the differences between natural magnetic fields and electromagnetic fields.

SCHULTZ's Law:

Weak stimuli have a sustaining effect on life processes,
 Moderate stimuli not only improve, but also exercise these processes,
 Strong stimuli have a negative effect, inhibiting reactions or causing paradoxical reactions.

12. FEELING OF CURRENT only on one electrode.

Some people may be concerned that they do not feel the current on both electrodes. This is due to the direction of the current flow and is also based on the laws of biophysics. It is normal and does not indicate any improper functioning of the device.

If the user really wants to feel the current on both electrodes and has experience, they can set the mode in the LIFE COMMANDER program to "BIPOLAR." Positive polarization is felt in the right hand, while negative polarization is felt in the left hand.

It is worth mentioning that the negative pole has a stimulating effect, while the positive pole has an analgesic effect. The negative pole is used in cases of sensory disorders or to prevent degeneration of nerve fibers in a damaged nerve. The positive pole has a negative impact on damaged nerves, which is why in such conditions, the cathode is always used as the active electrode. The effects of each electrode are presented below.

Positive polarity +	Negative polarity -
Causes a decrease in pH.	Causes an increase in pH (alkaline reaction).
Acts as a coagulant and thus causes tissue hardening.	Acts as a softening and dissolving agent.
Reduces tissue excitability.	Increases tissue excitability.
Causes a decrease in muscle tension.	Causes an increase in muscle tension. Acts vasodilatory (causes skin redness).
Affects slight dehydration of the skin.	Can cause slight skin swelling.
Acts soothing, toning, and calming.	Acts irritant.

CAUTION: "Feeling the electric current"

Sweeper is a device that operates based on resonant-frequency action. The main factor that influences its action is frequency, frequency sequence, and their temporal characteristics. Electric current or electromagnetic field is merely used for transmitting frequency information to the body.

We recommend setting the Intensity to a minimum. You can increase it if necessary, but only until you feel a tingling sensation, and then slightly decrease it. When working with induction, the field parameters are optimally adjusted by the manufacturer and do not require regulation.

It should be understood that "feeling the current" or the lack thereof is not an indicator of the device's proper functioning. These are purely individual and subjective feelings.

13. The Sweeper CT does not help to solve the problem for which it was purchased.

First, make sure that Sweeper CT is working. In many cases, simply checking the position of the electrodes and loops greatly improves the device's performance.

It is also possible that the duration of use was not long enough to achieve the desired effect. Consider using Sweeper CT regularly for 2-3 weeks. The beneficial effect is usually observed within the first few days.

It is important to know that we are always ready to communicate with you, and if you have any questions, please contact us through our social networks or email, also on our website.

14. The device has been exposed to excessive amounts of water.

The Sweeper CT was not designed to be waterproof, although short exposure to water should not permanently damage the device. If the Sweeper CT is exposed to water, recovery may be possible by immediately turning off the device and placing it in a sealed bag or container with uncooked dry rice for 2-3 days. Rice acts as a drying agent and slowly removes moisture.

15. The device does not turn on, the screen does not light up.

Charge the battery. Check the USB cable.

16. How to check the intensity of electromagnetic field induction.

Turn on the device. In the device menu, select "Settings", and then "Inductor Test". Place the field indicator (tester) provided with the device at the back of the device or near the external coil/inductive mat. If the LED on it lights up, it means that the built-in inductor or external loop/matt (if an external inductive loop is connected to the device, then the built-in inductor does not work at that time, and vice versa) is functioning properly. If you increase the distance, the brightness will gradually diminish. This happens because the magnetic field intensity decreases with distance from the surface of the inductor according to the inverse square law.

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